

Filtration System

Glass Solvent Vacuum Filters



Introduction

- ▶ This product is made of high-quality extra hard glass. It is sparking and crystal-clear. Its thickness is even.
- ▶ The performance of withstanding voltage and its leakproofness is very good. Its size agrees with the international standard size.

Order Information:

Item	Description	Packing
GVF30010	1) 300ml Glass Funnel with cover 2) 1000ml Glass Solvent Collection 3) Stainless Steel Clamp 4) Solid Glass Frit Packaging: 28cm*29cm*29cm	1set
GVF30020	1) 300ml Glass Funnel with cover 2) 2000ml Glass Solvent Collection 3) Stainless Steel Clamp 4) Solid Glass Frit Packaging: 32cm*33cm*24cm	1set
GVF50020	1) 500ml Glass Funnel 2) 2000ml Glass Solvent Collection 3) Stainless Steel Clamp 4) Solid Glass Frit Packaging: 32cm*33cm*24cm	1set

Filtration System

Multiple Vacuum Filtration System

Introduction

- ▶ Each station uses separate control valve for independent operation, easy to use and disinfect with a high efficiency.
- ▶ This system can filtrate three or six samples at the same time at low rate, with a high sensitivity.
- ▶ Sturdy won't tip when full loaded due to the units with low center of gravity.
- ▶ Anodized aluminium handles at both ends are designed for positioning on bench top.
- ▶ All MVFseries vacuum filtration devices are made of stainless steel which is suitable for filtration, especially in the microbiological field.

Part No.	Description	
MVF0101S	<ol style="list-style-type: none">1) Stainless Steel Manifold2) One 300ml Stainless Steel Filtration Units3) One Glass Collect Bottle4) Two pieces of hosepipe	
MVF0101G	<ol style="list-style-type: none">1) Stainless Steel Manifold2) One 300ml Glass Filtration Units3) One Glass Collect Bottle4) Two pieces of hosepipe	
MVF0103S	<ol style="list-style-type: none">1) Stainless Steel Manifold2) Three 300ml Stainless Steel Filtratio Units3) One Glass Collect Bottle4) Two pieces of hosepipe	
MVF0103G	<ol style="list-style-type: none">1) Stainless Steel Manifold2) Three 300ml Glass Filtration Units3) One Glass Collect Bottle4) Two pieces of hosepipe	
MVF0106S	<ol style="list-style-type: none">1) Stainless Steel Manifold2) Six 300ml Stainless Steel Filtration Units3) One Glass Collect Bottle4) Two pieces of hosepipe	
MVF0106G	<ol style="list-style-type: none">1) Stainless Steel Manifold2) Six 300ml Glass Filtration Units3) One Glass Collect Bottle4) Two pieces of hosepipe	

Part No.	Description	
MVFSM01	Stainless Steel Manifold, single	
MVFSM03	Stainless Steel Manifold, For Three Units	
MVFSM06	Stainless Steel Manifold, For Six Units	
MVFF300G	Glass Funnel, 300ml	
MVFF300S	Stainless Steel Funnel, 300ml	