


Filter Series

Reusable Syringe Filter Holder

Introduction

- ▶ Plastic syringe filter holder made of pure polypropylene.
- ▶ Stainless steel syringe filter holder made of 316L stainless steel.



Order Information

P/N	Description	Packing
SFH013P	Reusable Syringe Filter Holder, PP Dia. 13mm	10
SFH025P	Reusable Syringe Filter Holder, PP Dia. 25mm	10
SFH050P	Reusable Syringe Filter Holder, PP Dia. 50mm	10
SFH013S	Reusable Syringe Filter Holder, 316L Stainless Steel Dia. 13mm	1
SFH025S	Reusable Syringe Filter Holder, 316L Stainless Steel Dia. 25mm	1
SFH047S	Reusable Syringe Filter Holder, 316L Stainless Steel Dia. 47mm	1