


Liquid handling

Pipette Tip



Introduction

- ▶ Made of PP.
- ▶ Volume: 10ml, 200ml, 1000ml.
- ▶ Smooth finish for minimum retention of collected samples.
- ▶ Available in bulk package or packed in box.
- ▶ Non Sterile/Gamma Sterile/EO Sterile.

Order Information

Item No.	Description	Volume(μl)	Type	With Filter	Sterile	Package
MPT1010GN	Pipette Tips	10	Fit for Gilson	No	No	1000pcs/pk,50pks/cs
MPT1200GN	Pipette Tips	200	Fit for Gilson	No	No	1000pcs/pk,50pks/cs
MPT11000GN	Pipette Tips	1000	Fit for Gilson	No	No	500pcs/pk,30pks/cs
MPT1010EN	Pipette Tips	10	Fit for Eppendorf	No	No	1000pcs/pk,50pks/cs
MPT1200EN	Pipette Tips	200	Fit for Eppendorf	No	No	1000pcs/pk,50pks/cs
MPT11000EN	Pipette Tips	1000	Fit for Eppendorf	No	No	500pcs/pk,30pks/cs
MPT0010GN	Pipette Tips	10	Fit for Gilson	Yes	No	1000pcs/pk,50pks/cs
MPT0200GN	Pipette Tips	200	Fit for Gilson	Yes	No	1000pcs/pk,50pks/cs
MPT01000GN	Pipette Tips	1000	Fit for Gilson	Yes	No	500pcs/pk,30pks/cs
MPT2010N	Pipette Tips	10	Fit for Gilson	No	No	96pcs/boxes, 120boxes/cs
MPT2200N	Pipette Tips	200	Fit for Gilson	No	No	96pcs/boxes, 100boxes/cs
MPT21000N	Pipette Tips	1000	Fit for Gilson	No	No	100pcs/box60boxes/cs
MPT3010N	Pipette Tips	10	Fit for Gilson	Yes	No	96pcs/boxes, 120boxes/cs
MPT3200N	Pipette Tips	200	Fit for Gilson	Yes	No	96pcs/boxes, 100boxes/cs
MPT31000N	Pipette Tips	1000	Fit for Gilson	Yes	No	100pcs/box60boxes/cs
MPT2010S	Pipette Tips	10	Fit for Gilson	No	Yes	96pcs/boxes, 120boxes/cs
MPT2200S	Pipette Tips	200	Fit for Gilson	No	Yes	96pcs/boxes, 100boxes/cs
MPT21000S	Pipette Tips	1000	Fit for Gilson	No	Yes	100pcs/box60boxes/cs
MPT3010S	Pipette Tips	10	Fit for Gilson	Yes	Yes	96pcs/boxes, 120boxes/cs
MPT3200S	Pipette Tips	200	Fit for Gilson	Yes	Yes	96pcs/boxes, 100boxes/cs
MPT31000S	Pipette Tips	1000	Fit for Gilson	Yes	Yes	100pcs/box60boxes/cs

PS: Gamma sterife is available.



Liquid handling

Pasteur Pipette



Introduction

- ▶ Made of LDPE.
- ▶ Volume: 1ml, 3ml, 5ml, 10ml.
- ▶ Suitable for liquid transfer.

Order Information

Item No.	Description	Size	Material	Package
MPP001	Pasteur Pipette	1ml	LDPE	5000pcs/cs
MPP003	Pasteur Pipette	3ml	LDPE	5000pcs/cs
MPP005	Pasteur Pipette	5ml	LDPE	5000pcs/cs
MPP010	Pasteur Pipette	10ml	LDPE	5000pcs/cs

Serological Pipette



Introduction

- ▶ Made of high-grade polystyrene, excellent for clear observation and reduce liquid attachment on the pipet.
- ▶ Volume: 1ml, 2ml, 5ml, 10ml, 25ml, 50ml. Graduations Inaccuracy < 2% at full volume.
- ▶ Color-coded rings and package for easier identifying the correct size.
- ▶ All pipets are supplied with a filter plug to prevent overflow.
- ▶ Available with gamma sterile and non sterile package.
- ▶ Non-pyrogenic, Non-cytotoxic and Non-hemolytic.

Order Information

Item No.	Size(ml)	Graduation	Package	Gamma Sterile	PCS/CTN
MLSP001S	1	1/100ml	Individually pack, 50pcs/bag, 10bags/box	Yes	3000
MLSP002S	2	1/100ml	Individually pack, 50pcs/bag, 10bags/box	Yes	3000
MLSP005S	5	1/10ml	Individually pack, 50pcs/bag, 10bags/box	Yes	1500
MLSP010S	10	1/10ml	Individually pack, 50pcs/bag, 8bags/box	Yes	1200
MLSP025S	25	2/10ml	Individually pack, 25pcs/bag, 8bags/box	Yes	600
MLSP050S	50	5/10ml	Individually pack, 25pcs/bag, 4bags/box	Yes	300
MLSP001N	1	1/100ml	50pcs/bag, 10bags/box	No	3000
MLSP002N	2	1/100ml	50pcs/bag, 10bags/box	No	3000
MLSP005N	5	1/10ml	50pcs/bag, 10bags/box	No	1500
MLSP010N	10	1/10ml	50pcs/bag, 8bags/box	No	1200
MLSP025N	25	2/10ml	25pcs/bag, 8bags/box	No	600
MLSP050N	50	5/10ml	25pcs/bag, 4bags/box	No	300

Tubes

Centrifuge Tube



Introduction

- ▶ Made of high quality, crystal-clear polypropylene, chloroform-resistant, leak-proof.
- ▶ Widely used in chemical and life science research, suitable for most brand of centrifuge.
- ▶ Available with almost all kinds of volume, with conical bottom and self-standing bottom, with white printing area for easy mark.
- ▶ Autoclavable at 121°C and freezable to -80°C.

Order Information

Item No.	Volume(ml)	Type	Sterile	Package
CT000NB	0.2	Conical	No	1000pcs/pk,70pks/cs
CT000NA	0.5	Conical	No	1000pcs/pk,30pks/cs
CT001NA	1.5	Conical	No	500pcs/pk,20pks/cs
CT002NA	2	Conical	No	500pcs/pk,20pks/cs
CT005NA	5	Round	No	500pcs/pk,20pks/cs
CT007NA	7	Round	No	300pcs/pk,20pks/cs
CT010NA	10	Round	No	100pcs/pk,30pks/cs
CT010NB	10	Round	No	100pcs/pk,30pks/cs
CT015NA	15	Conical	No	100pcs/pk,5pks/cs
CT015SA	15	Conical	Yes	100pcs/pk,5pks/cs
CT050NA	50	Conical	No	50pcs/pk,10pks/cs
CT050SA	50	Conical	Yes	50pcs/pk,10pks/cs
CT150NA	50	Self-standing	No	50pcs/pk,10pks/cs
CT150SA	50	Self-standing	Yes	50pcs/pk,10pks/cs

Tubes

Test Tube



Introduction

- ▶ Made of PP and PS.
- ▶ Available with round bottom and conical bottom.

Order Information

Item No.	Description	Size(mm)	Material	Package
TTR1260S	Round Bottom	12*60	PS	6000pcs/cs
TTR1260P	Round Bottom	12*60	PP	6000pcs/cs
TTR1275S	Round Bottom	12*75	PS	5000pcs/cs
TTR1275P	Round Bottom	12*75	PP	5000pcs/cs
TTR1375S	Round Bottom	13*75	PS	5000pcs/cs
TTR1375P	Round Bottom	13*75	PP	5000pcs/cs
TTR1378S	Round Bottom	13*78	PS	5000pcs/cs
TTR1378P	Round Bottom	13*78	PP	5000pcs/cs
TTR13100S	Round Bottom	13*100	PS	4000pcs/cs
TTR13100P	Round Bottom	13*100	PP	4000pcs/cs
TTR15100S	Round Bottom	15*100	PS	3000pcs/cs
TTR15100P	Round Bottom	15*100	PP	3000pcs/cs
TTR16100S	Round Bottom	16*100	PS	2500pcs/cs
TTR16100P	Round Bottom	16*100	PP	2500pcs/cs
TTC16102S	Conical Bottom	16*102	PS	2500pcs/cs
TTC16102P	Conical Bottom	16*102	PP	2500pcs/cs
CTT1275S	Cell Test Tube	12*75	PS	5000pcs/cs
CTT15100S	Cell Test Tube	15*100	PS	5000pcs/cs
TTR22120S	Round Bottom	22*120	PS	1200pcs/cs
TTR22153S	Round Bottom	20*153	PS	1200pcs/cs
TS12	Tube Stopper	φ 12	LDPE	30000pcs/cs
TS13	Tube Stopper	φ 13	LDPE	30000pcs/cs
TS15	Tube Stopper	φ 15	LDPE	15000pcs/cs
TS16	Tube Stopper	φ 16	LDPE	15000pcs/cs